

D. CLONE-FISH

Here we describe a method using clones for adjusting the hybridization conditions. In brief, clones carrying the target sequence of the new probe are grown with chloramphenicol and IPTG. If a vector with an inducible promoter upstream of the multiple cloning site is chosen then this treatment leads to an *in vivo* transcription of the cloned 16S rDNA and accumulation of 16S rRNA of the uncultured organism inside the *E. coli* cell. After standard fixation of these clones they could be used as analogs to cultured organisms for determining the melting point of probes.

please cite for this protocol: Schramm, A., B. M. Fuchs, J. L. Nielsen, M. Tonolla, and D. A. Stahl. 2002. Fluorescence *in situ* hybridization of 16S rRNA gene clones (Clone-FISH) for probe validation and screening of clone libraries. *Environmental Microbiology* 4:713-720.

Theory

Chloramphenicol:

- blocks protein synthesis at the ribosome;
- Rop -protein is not produced any more => leads to uncontrolled (increased) plasmid copy number (for ColE1 origin plasmids like pMB1, pGEMt, pBR322, pUC, pBLUESCRIPT)
- cell elongation

IPTG:

- induces transcription of T7-polymerase => in turn transcribes the insert from T7 promoter on plasmid
- induces transcription of lacZ gene => insert in opposite direction than T7 transcribed
- T7 has higher efficiency than lacZ

Protocol

1. select appropriate vector-host combination [e.g. pGEM-T in JM109(DE3) (for easy A-T-cloning of PCR products) or pET-23(+) in NovaBlue(DE3) (for blunt- or sticky-end cloning)];
2. grow the clones to OD600 of 0.3-0.4 from 1:1000 dilution of an overnight culture;
3. add IPTG (1 mM) and incubate for 1 h;
4. add chloramphenicol (170 mg l⁻¹) and incubate for 4 h;
5. fix with paraformaldehyde and store in PBS:ethanol (1:1) at -20 °C

Clones can be used for FISH for more than one year.

Note: only probes that do not target rRNA of the *E. coli* host cells are suitable for Clone-FISH.